

Diversey

Oxivir® Plus
Disinfectant
Cleaner
Concentrate



Dilutable one-step disinfectant cleaner, hydrotherapy tub cleaner and sanitizer, formulated with proprietary hydrogen peroxide technology. Now Ecologo certified.

Virucidal and Bactericidal

- Disinfects in just five minutes
- Effective against Human Immunodeficiency Virus (HIV-1) (AIDS Virus), Polio virus, Pseudomonas aeruginosa, Staphylococcus aureus, Methicillin resistant Staphylococcus aureus (MRSA), Salmonella choleraesuis, Influenza A Virus and VRE
- Kills Norovirus and Canine Parvovirus
- Bactericidal 5 minutes
- Tuberculocidal 5 minutes
- Virucidal 5 minutes
- Fungicidal 5 minutes
- Non-Food Contact Sanitizer 30 seconds at 1:128

Compatible With Hard Surfaces

- Active breaks down into water and oxygen by-products
- No rinsing required
- Delivers fast and effective cleaning
- Not irritating to eyes and skin in use dilution

Recommended For

- Healthcare
- Education
- Retail
- BSC
- Government







Use Overview



Oxivir® Plus Disinfectant Cleaner Concentrate is a one-step disinfectant cleaner using patented AHP® technology for use on walls, stainless steel, glazed porcelain, plastic surfaces, hydrotherapy tubs and other hard, nonporous surfaces.

- To disinfect, apply a 1:40 dilution to thoroughly wet surface.
- Let solution remain on surface for a minimum of 5 minutes.
- Wipe or allow to air dry.
- Prepare a fresh solution if the solution becomes visibly dirty or diluted.

Dilutable Disinfectant Cleaner

- 1:40 dilution Bactericidal and virucidal, including antibiotic-resisting bacteria, in 5 minutes.
- 1:128 dilution Non-food contact sanitizer in 30 seconds or general purpose cleaning.



Product Specifications

Description	Concentrate
Color/Form	Colorless liquid
DIN Reg. No.	02403684
Flash Point	> 200° F (> 93.4° C)
-11	0.8 (Concentrate)
pH 	1.9 (Use Dilution 1:40)
Scent	Characteristic
Shelf Life	2 years
	30 days (Diluted)
Solubility	Complete
Specific Gravity	1.037

Hazard Rating

HMIS	Concentrate	NFPA
2	Health	2
0	Flammability	0
0	Reactivity	0
HMIS	Use Dilution 1:16	NFPA
HMIS 0	Use Dilution 1:16 Health	NFPA 0
0 0		NFPA 0 0

4=Very High; 3=High; 2=Moderate; 1=Slight; 0=Insignificant.

Available Items

Product Code	Description/Package Size	Dilution
5919024	4 x 1 gallon / 3.78 L Containers	1:40, 1:128
5919032	2 x 1.5 L RTD® Units	1:40
5919059	2 x 84.5 oz. / 2.5 L J-Fills®	1:40

Safety Reminder

Please make sure your employees read and understand the product label and Material Safety Data Sheet before using this product. The label contains directions for use, and both the label and MSDS contain hazard warnings, precautionary statements and first aid procedures. MSDS are available on-line at www.diversey.com or by calling 888 352 2249.

Reference Sheet Fiche de référence



Oxivir® Plus Disinfectant Cleaner Concentrate Nettoyant Désinfectant concentré Oxivir® Plus



For Use in Animal Housing Areas, Medical and Healthcare Establishments

A POWERFUL DISINFECTANT EFFECTIVE AGAINST:

BACTERICIDAL:

Staphylococcus aureus (ATCC 6538) Pseudomonas aeruginosa (ATCC 15442) Salmonella enterica (ATCC 10708) Klebsiella pneumonià (ATCC 13882) Enterobacter aerogenès (ATCC BAA-2356) Methicillin-resistant *Staphylococcus aureus* (MRSA) (ATCC 29247) Vancomycin-resistant *Enterococcus faecalis* (VRE) (ATCC 51299) Escherichia coli 0157:H7 (ATCC 43865)

TUBERCULOCIDAL:

Mycobacterium terrae (ATCC 15755)

FUNGICIDAL:

Trychophyton mentagrophytes (ATCC 9533) Candida albicans (ATCC 10231)

VIRUCIDAL:

Proven effectiveness against the Poliovirus Type 1, Sabin strain type 1 (ATCC VR-192) which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses. HIV-1 Human Immunodeficiency virus (HIV), Strain HTLV-IIIB (HIV-1)
Feline Calicivirus, F9 Strain (ATCC VR-782), as a surrogate for Norwalk and Norwalk-like viruses
Rhinovirus serotype 14 (ATCC VR-284)
Canine Parvovirus (CPV), the Cornell strain (ATCC VR-2017)

Influenza Virus PR8 strain (ATCC VR-95)
Rotavirus Wa strain (ATCC VR-2018)
Murine Norovirus type 1 (MNV-1) strain S99

Human Coronavirus 229E (ATCC VR-740)

Human Parainfluenza Virus 3 (ATCC VR-92)
This product has demonstrated effectiveness against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus.

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZING:

Staphylococcus aureus (ATCC 6538) Pseudomonas aeruginosa (ATCC 15442) Salmonella enterica (ATCC 10708) Klebsiella pneumonia (ATCC 13882 Campylobacter jejuni (ATCC 33560

Methicillin-resistant *Staphylococcus aureus* (MRSA) (ATCC 29247) Vancomycin-resistant *Enterococcus faecalis* (VRE) (ATCC 51299) Acinetobacter baumanii (ATCC 19606) Escherichia coli 0157:H7 (ATCC 43865)

VIRUCIDAL: 5 min. at 20° C TUBERCULOCIDAL: 5 min. at 20° C FUNGICIDAL: 5 min. at 20° C

BROAD-SPECTRUM SANITIZER: 30 sec. at 20° C

See container label for First Aid, Precautionary Statements and Directions for Use

Pour une utilisation dans les zones d'hébergement des animaux, les établissements médicaux et de santé

UN PUISSANT DÉSINFECTANT EFFICACE CONTRE :

BACTÉRICIDE:

Staphylococcus aureus (ATCC 6538) Pseudomonas aeruginosa (ATCC 15442) Salmonella enterica (ATCC 10708) Klebsiella pneumonia (ATCC 13882) Enterobacter aerogenes (ATCC BAA-2356)

Staphylococcus aureus résistant à la méthicilline (MRSA) (ATCC 29247) Enterococcus faecalis résistant à la vancomycine (VRE) (ATCC 51299) Acinetobacter baumanii (ATCC 19606) Escherichia coli 0157:H7 (ATCC 43865)

TUBERCULOCIDE:

Mycobacterium terrae (ATCC 15755)

FONGICIDE:

Trychophyton mentagrophytes (ATCC 9533), Candida albicans (ATCC 10231)

Efficacité prouvée contre le Poliovirus Type 1, Sabin strain type 1 (ATCC VR-192) ce qui permet une réclamation à large spectre contre la plupart virucide virus enveloppés et non enveloppés viruses.

Virus VIH-1 de l'immunidéficience humaine (VIH), souche HTLV-IIIB (HIV-1), calicivirus félin, souche F9 (ATCC VR-782), comme un substitut pour les rhinovirus Norwalk et similaires sérotype 14 (ATCC VR-284), parvovirus canin(CPV), souche Cornell (ATCC VR-2017), virus grippal souche PR8 (ATCC VR-95), Souche de rotavirus Wa (ATCC VR-2018), Norovirus murin de type 1 (NVM-1) souche S99, Coronavirus humain 229E (ATCC VR-740), Virus parainfluenza humain 3 (ATCC VR-92)

Ce produit a prouvé son efficacité contre le Poliovirus et est prévu de désactiver tous les virus de la grippe A dont le virus grippal A pandémique 2009 (H1N1).

ASSAINISSANT À LARGE SPECTRE SUR DES SURFACES ENVIRONNEMENTALES SANS CONTACT AVEC DES ALIMENTS :

Staphylococcus aureus (ATCC 6538) Pseudomonas aeruginosa (ATCC 15442) Salmonella enterica (ATCC 10708) Klebsiella pneumonia (ATCC 13882) Campylobacter jejuni (ATCC 33560)

Staphylococcus aureus résistant à la méthicilline (MRSA) (ATCC 29247) Enterococcus faecalis résistant à la vancomycine (VRE) (ATCC 51299) Acinetobacter baumanii (ATCC 19606) Escherichia coli 0157:H7 (ATCC 43865)

VIRUCIDE: 5 Min. à 20 °C TUBERCULOCIDE: 5 Min. à 20 °C FONGICIDE: 5 Min. à 20 °C BACTÉRICIDE: 5 Min. à 20 °C

DÉSINFECTANT À LARGE SPECTRE D'EFFICACITÉ : 30 sec. à 20 °C

Voir l'étiquette apposée sur le récipient pour les premiers secours, les énoncés de précaution et les instructions d'utilisation

Manufactured for / Fabriqué pour : ©2015 Diversey, Inc., 8310 16th Street, Sturtevant, WI 53177-1964 Distributed by / Distribué par : Diversey Canada, Inc., 3755 Laird Road, Units 8 to 11, Mississauga, ON, Canada L5L 0B3

REF90223 (15/041)

NOTE: Read and check carefully all information, including the specifications in this box. The job will not be sent to print until final okay is received, All information contained on the label and in this box are the responsibility of the final approver.



Black

Brand: Sealed Air / Diversey

Product: Oxivir Plus

Pack Size: Reference Sheet

Date Started: 02/06/15

Size: 8.5" x 11"

Product Code: 5919024

Tracking: REF90223 Replaces: REF90167

Minimum type: 6/6

Bar Code Type: UPC-A

Coordinator: Jennifer Keszler

Designer: D. Reeves

Revision History: 02/06/15: Update/dr 02/10/15: Bolded Text/dr