

# Sporicidal Hard Surface Disinfectant



---

Ready-to-use liquid sporicide for hard surfaces that is effective against *Clostridium difficile*, *Clostridium sporogenes* and *Bacillus subtilis*, utilizing AHP® Hydrogen Peroxide technology. Hard surface disinfectant against polio, pseudomonas, staph and salmonella; fungicidal against trichophyton.

## Effective

- ▶ Ready-to-use liquid sporicide for hard surfaces
- ▶ Effective against *Clostridium difficile*
- ▶ No mixing needed
- ▶ 2 year shelf life

## Easy-to-Use

- ▶ Just apply and wipe

---

## Designed For

- ▶ Healthcare
- ▶ Hospitality



# Use Overview

**Sporicidal Hard Surface Disinfectant** is a ready-to-use liquid sporicide for hard surfaces that is effective against *Clostridium difficile*, *Clostridium sporogenes* and *Bacillus subtilis*, utilizing AHP® Hydrogen Peroxide technology. Hard surface disinfectant against polio, pseudomonas, staph and salmonella; fungicidal against trichophyton.

## Sporicidal Disinfection of Hard, Nonporous Surfaces Use Directions

- ▶ This product is not to be used as a terminal Sterilant/High Level Disinfectant for critical instruments and medical devices. This product is intended for use as an Environmental Surface Sporicide on hard, nonporous surfaces in C. difficile patient rooms. (See Notice to User.)
- ▶ Apply liquid to surfaces using a cloth or disposable wipe. To achieve Sporicidal disinfection, ensure surface remains wet for 10 minutes at 20°C.
- ▶ Thoroughly rinse to ensure all product residue has been removed.

## Personal Protection

- ▶ Disposable gloves, gowns, face masks, or eye coverings as appropriate, must be worn during all cleaning of body fluids, and decontamination procedures.

## Cleaning Surfaces Contaminated With Body Fluids

- ▶ Appropriate personal protection must be worn during all cleaning and disinfection procedures. Remove excess body fluids with absorbent materials.
- ▶ Apply liquid to surfaces using a cloth or disposable wipe and let stand for 1 minute, thoroughly rinse.
- ▶ Disinfect contaminated area.

**NOTICE TO USER: Sporicidal Hard Surface Disinfectant** is not intended for use on medical equipment without device manufacturer approval.

## Disposal of Infectious Material

- ▶ Products contaminated with blood or body fluids should be disposed of according to Federal, Provincial and local regulations for infectious waste disposal.



## Product Specifications

Description	
<b>Color/Form</b>	Clear colorless liquid
<b>DIN</b>	02376008
<b>Flash Point</b>	>200° F (>93.4° C)
<b>pH</b>	2.7
<b>Scent</b>	Characteristic
<b>Shelf Life</b>	2 years
<b>Solubility</b>	Complete
<b>Specific Gravity</b>	1.0325

## Hazard Rating

HMIS	Concentrate	NFPA
2	Health	2
0	Flammability	0
0	Reactivity	0

4=Very High; 3=High; 2=Moderate; 1=Slight; 0=Insignificant.

## Available Items

Product Code	Description/Package Size	Dilution
5728148	12 x 32 oz./946 mL Squeeze bottles	Ready-to-Use

## Safety Reminder

Please make sure your employees read and understand the product label and Material Safety Data Sheet before using this product. The label contains directions for use, and both the label and MSDS contain hazard warnings, precautionary statements and first aid procedures. MSDS are available on-line at [www.diversey.com](http://www.diversey.com) or by calling 888 352 2249.

AHP® Technology and Logos are registered trademarks of Virox Technologies Inc., used under license.

For additional information, please contact us at +1 800 558 2332 or at [www.diversey.com](http://www.diversey.com)

© 2012 Diversey, Inc., 8310 16th Street, Sturtevant, WI 53177-1964. All Rights Reserved. SPC90175 (12/115)

