

Product: CEDERROTH WALL BRACKET FOR EYE WASH
Reference code: 7200



1. General description

A holder in green hard plastics intended for storage of a Cederroth Eye Wash bottle, making it easy to find in an emergency situation. The holder can be mounted on any vertical surface e.g. on a wall, in a vehicle, at the laboratory bench or on the Cederroth First Aid Panel. An afterglowing sign with instructions for use is included.

The Cederroth Wall Bracket for eye wash also serves as an **opening device** for the Cederroth Eye Wash and as an **anti-theft system**, as the Eye Wash bottle is automatically opened when removed from the holder and thus can not be returned into the holder.



2. Intended use

The Wall Bracket is intended for providing a safe and organized storage of the Eye Wash products, which will result in easy access of the products whenever needed.

3. Technical description

Product dimensions :

Wall bracket (l x w x h) mm: 84 x 262 x 70
 Sign (l x w x h) mm: 85 x 250

Material

Wall bracket Polyamid plus Glassfibre
 Afterglowing sign PVC

4. Composition

Consumer package

Transport package

| <u>Ref#</u> | <u>Product/Material</u> | <u>Pcs</u> | Number of consumer packages | |
|-------------|-------------------------|------------|-----------------------------|-----------------|
| 51616100 | Wall bracket | 1 | in transport package: | 10 |
| 51672170 | Afterglowing sign | 1 | Dimensions (l x w x h) mm: | 390 x 290 x 200 |
| 51618671 | Plastic zip bag | 1 | Weight (kg): | 2,3 |



PRODUCT SPECIFICATION

5. Labeling

On product/ consumer package

Cederroth address, telephone- and faxnumber, web-address

Product name, Reference code

Instruction for use

EAN-code

On Wall Bracket:

Lot-number

Reference code for Eye Wash refill

Symbols: LOT

On transport package

Cederroth address, telephone- and faxnumber, web-address

Product name, Reference code

Lot-number, EAN-code

Quantity

Symbols: LOT

6. Storage conditions

The product should be stored in the original package for maximum protection until time of use.

7. Accessories

The wall bracket can be mounted on the Cederroth First aid panel.

Authorization: